


New and Improved COVID-19 Self-Collection Test


Nasal collection – front of nose



- Position head slightly back, and insert the swab into the LEFT nostril, only about one inch.
- Rotate the swab in a circular motion around the entire inside edge of the nostril. Do this 2 times and then keep it in place for 15 seconds.
- Repeat the same process in the RIGHT nostril using the same swab.

Nueva y Mejorada Prueba de Autohisopado de COVID-19

Toma de muestras nasales — frente de la nariz



- Incline la cabeza ligeramente hacia atrás e inserte el hisopo en la fosa nasal **IZQUIERDA**, solo una pulgada (2.5 cm) aproximadamente.
- Gire el hisopo con un movimiento circular por todo el borde interior de la fosa nasal. Hágalo 2 veces y déjelo allí durante 15 segundos. Repita el mismo proceso en la fosa nasal **DERECHA** con el mismo hisopo.